




D2D ARTWORK PREPRESS MASTERFILE	Product Brand	ThickenUp	Selling Markets	DE-AT-CZ-SK
	Denominator	Clear	Production Plant	Osthofen - DE
	Packaging size / weight	125g	Tech. drawing / spec nbr.	—
	Packaging type	Tin-Paper	Printer - Print process	Druckerei Plenck KG - TBC
	Variant	composite	Colour Profile	FOGRA 39L Mars
	ZPCK	44339684	Substrate color	Paper glued on composite can
	GTIN	7613032719616	Cap information	White
	FERT	12114002	Fonts	DIN Pro
	TD Code	516016	Min Font Size	6.75 pts – Height: 1.2 mm (E.U. Ready)
	Recipe ID	MD009	D2D Execution Manager	Aude Nicolino

File Name: ThickenUp_Clear_Neutral_Tin-Paper_125g_DE-AT-CZ-SK_AWP_44339684_TD516016

